precision®

moisture-displacing for highly-detailed impressions

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Illustrations for Jumbo Cartridge

Precision® Product Description

- Precision Extra Lite Body (Orange) ISO 4823, Type 3: Low Viscosity
- Precision Lite Body (Green) ISO 4823, Type 3: Low Viscosity
- Precision Medium Body (Blue) ISO 4823, Type 2: Medium Viscosity
- Precision Heavy Body (Purple) ISO 4823, Type 1: Heavy Viscosity
- Precision Putty (Purple) ISO 4823, Type 0: Putty Viscosity

Precision Technical Specifications

- Cartridge Mixing Technique (XL, L): Use a 1:1/1:2 dispensing cartridge gun, with static mixing tip provided (4.2-12D Diameter).
- Cartridge Mixing Technique (M, H): Use a 1:1/1:2 dispensing cartridge gun, with static mixing tip provided (6.5-11D Diameter), or 5:1 automix unit and blue dynamic mixing tip for Jumbo Cartridges.

Precision Extra Lite Body Color: Orange

- Linear dimensional change: max. -0.1 %
- Recovery from deformation: min. 99.6 %
- Strain in compression: <4.5 %
- Shore A Hardness: 44-49

Precision Lite Body Color: Green

- Linear dimensional change: max. -0.1 %
- Recovery from deformation: min. 99.6 %
- Strain in compression: <4.0 %
- Shore A Hardness: 46-51

Precision Medium Body Color: Blue

- Linear dimensional change:
- Recovery from deformation:
- Strain in compression:
- Shore A Hardness:

Cartridge
max. -0.2 %
min. 99.6 %
<3.5 %
58-62

Jumbo Cartridge
max. -0.25%
min. 99.6%
<4%
56-62

Precision Heavy Body Color: Purple

- Linear dimensional change:
- Recovery from deformation:
- Strain in compression:
- Shore A Hardness:

Cartridge

max. -0.2%
min. 99.3 %
max. 3.5%
59-63

Jumbo Cartridge

max. -0.25%
min. 99.3%
max. 3.5%
57-63

Precision Putty & Putty-Pak® Color: Purple

- Linear dimensional change: max. -0.2%
- Recovery from deformation: min. 99.0%
- Strain in compression: <2.5%
- Shore A Hardness: 68-72

- Storage conditions: In a dry place, at 18 –25 °C (64 –77 °F) and relative humidity of 50% or less. Do not expose to direct heat or sun. Store cartridges upright.
- Electroplating: Impressions may be copper or silver plated.
- Casting: The impression can be poured after 1 hour.

Indications for Use

Precision is an addition-reaction silicon impression material that is ideal for dual viscosity impression techniques, both one-step and two-step, as well as for a variety of dental procedures including, but not limited to, crowns, bridges, inlays, partial and complete dentures, denture repairs, implants and for forming functional peripheries.


Working and Setting Times

Precision offers a variable working time and setting reaction that is temperature controlled. If placed in the mouth immediately after extrusion, Precision will allow 45 seconds of intra-oral working time. If more time is required, Precision will allow for 1:45 extra-oral working time. Therefore, depending on how much working time is required (45 seconds – 1:45), the total set time will range from 2:45 – 3:45 minutes.

Directions for Use


Instructions for cartridge:

1. Insert the cartridge into the cartridge dispensing gun.
2. Remove the cartridge closure.
3. Before applying the mixing tip, extrude impression material until it exits both chambers at the same rate.
4. Attach the mixing tip on the cartridge and extrude the required amount of material directly onto the prepared tooth or into an impression tray or delivery syringe.
5. Leave the filled mixing tip on the cartridge after use.
6. Before using the cartridge again, remove and discard the used mixing tip. Before re-applying another mixing tip, be sure to extrude material so the product is exiting both chambers at the same rate.
7. Attach a fresh mixing tip and proceed as usual.


Tray Filling Time:

Up to 1:45 min.
(73.4° F / 23° C)

Intraoral Syringe Time:

Up to 45 sec.
(95° F / 35° C)

Total Set Time:

2:45 – 3:45 min.
(73.4° F / 23° C)

Instructions for Jumbo Cartridge

Please refer to the illustrations in the front of this booklet.

1. Prior to initial use of product, remove the blue stopper cap by turning it slightly to the left while pulling it out simultaneously. Discard cap after removal (fig. 1).
2. Insert cartridge into the automix dispensing machine.
3. Attach the blue dynamic mixing tip to the end of the cartridge (fig. 2). To lock tip in place, line up outlets, insert and turn to the right.
4. Start machine and extrude about $\frac{1}{4}$ inch of material to ensure a homogenous mix of base and catalyst. Discard this material (fig. 3).
5. Fill the tray with the desired amount of material by holding the tray in a tilted position and gently pressing against the dynamic mixer. Keep the tip immersed in the material to avoid air bubbles (fig. 4).
6. Do not remove the used mixing tip from the cartridge body as it acts as a seal.
7. Before using again, remove used tip and check that the outlets in the activating heads are not blocked. Fit a new mixing tip in place and continue with steps 4 – 6.

Instructions for Putty Jars

1. Scoop even amounts of base and catalyst from jars.
2. Make sure to keep base and catalyst scoops separate from each other.
3. Knead putty thoroughly until the color is homogenous and streak free.
Note: Do not mix while wearing latex gloves; vinyl or nitrile gloves are recommended.
4. Load mixed putty into tray.

Instructions for Putty-Pak®

1. Peel off lid from Putty-Pak.
2. Fold base and catalyst together.
3. Squeeze product together.
4. Remove product from Pak and mix thoroughly until color is homogenous and streak free.
Note: Do not mix while wearing latex gloves; vinyl or nitrile gloves are recommended.
5. Load mixed putty into tray.

Special Notes & Precautions

- Cartridge material must be used with a 1:1/1:2 dispensing cartridge gun system. Jumbo Cartridge material can be used in a 5:1 compatible automix unit with dynamic mixer tips.
- Solutions applied prior to impression taking (such as retraction solution) may impede the setting reaction of the impression. Therefore, material must be cleaned off thoroughly with water and dried.
- Do not leave remains of the impression material in the sulcus.
- Avoid splattering this material onto clothing as it cannot be removed by dry cleaning.
- The possibility of allergic reactions to the product in sensitive persons cannot be fully excluded. Possible reactions may include temporary numbness or tingling in the soft tissue.
- When taking an impression adjacent to an existing fixed bridge, make sure all undercuts are sufficiently blocked out.
- All intraoral areas touched by latex gloves should be rinsed and dried prior to impression taking.

The information provided by DenMat is based on comprehensive research and experience in application technology. Results are furnished to the best of our knowledge, subject to technical changes within the framework of product development. However, users must comply with and consider all recommendations and information in connection with any use.

Compatibility

Precision is an addition-reaction type silicone impression material and should not be combined with standard or condensation type silicon materials or polyether materials.

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Den-Mat Holdings, LLC
MDSS GmbH