Instruction Guidelines

Warnings and Precautions

DANGER: Not for use in the presence of flammable anesthetics.

WARNING: The Vitality Scanner 2006 should NOT be used on a patient with a pacemaker.

WARNING: Use of other accessories which are not authorized for use in connection with the Vitality Scanner 2006 may cause malfunction and may compromise patient safety.

CAUTION: This device has been investigated with regard to safety from electrical shock and fire hazard as well as electromagnetic compatibility (EMC). The device has not been investigated for other physiological effects.

CAUTION: For use by qualified and trained personnel only.

Package Contents

P/N	Description
973-0234	Vitality Scanner 2006 One (1) Short Tip One (1) Long Tip Grounding Lead with Lip Clip

Available Accessories

973-0232	Short Tip (3-pack)
973-0233	Long Tip (3-pack)
973-0231	Mini Crown Test Tip (2-pack)
973-0226	Grounding Lead with Lip Clip
973-0227	Lip Clip (5)

Symbols

Bottom Panel Symbols

Attention! Consult the manual before use.

CSA mark with "C/US" indicator for certified products.

MC 157498
(LR 113060)

CSA file number.

Type BF Patient Applied Part

Packaging / Handling Symbols

This way up.

Fragile, handle with care.

Keep away from rain.

Temperature limits for storage and transportation.

Index: page

English	1
Deutsch	3
Français	6
Italiano	8
Dansk	11
Español	13
Suomi	16
Nederlands	18
Svenska	21
Português	23
Norsk	26
Ελληνικά	28
日本	31

Product Features

The features of the SybronEndo Vitality Scanner 2006 are presented below in the form of a series of statements describing the function of the instrument in a step by step manner. Each statement is followed by a paragraph with related information.

1. When the probe makes contact with a tooth, the Vitality Scanner turns on automatically. The Vitality Scanner will turn on automatically as soon as you establish a good electrical contact between the probe tip and a tooth. An indicator light on the probe wand lights up when the Vitality Scanner turns on. The automatic turn-on feature prevents you from getting a false diagnosis due to poor electrical contact and it eliminates the possibility of jolting your patient by suddenly establishing good contact after initially having had poor contact between the probe and the tooth. The Vitality Scanner will not turn on when you test a tooth with a porcelain crown—a situation in which electrical pulp testing in general becomes invalid. An optional Mini-Tip is available which allows you access to the tooth below the crown margin. You may also test metallic fillings. Be cautious, however, with the readings you may receive from metallic fillings since the resistance of a metallic filling is considerably lower than that of natural tooth material, resulting in lower response levels.

2. After the unit turns on, the intensity of the electrical stimulus will increase automatically. The rate of increase in the strength of the stimulus can be controlled by a panel wheel. The voltage output is generated electronically so that the level rises slowly in the beginning of the test and more quickly at the end.

3. When your patient senses the stimulus, lift the probe from the tooth and read the stimulus level on the display. The patient will perceive the stimulation as a pulsating pressure, warmth, or tingling which slowly builds in strength unlike the sudden jolt common with other pulp testers. The display will stop incrementing when you break the contact between the probe and the tooth. The output voltage will reset to its lowest value after a delay of two seconds, but the display will continue to show the reading until you begin testing another tooth. If you reestablish contact within the delay interval, the pulp testing sequence will continue without resetting. This way you can avoid accidentally resetting the Vitality Scanner if you momentarily lose contact between the probe and the tooth.

4. The display will automatically reset when you move to another tooth. The automatic reset feature allows you to quickly and easily test a multitude of teeth simply by moving the probe from one tooth to the next. The output stimulus level resets automatically after each test. Unlike other pulp testers, our unit's reset feature prevents you from accidentally shocking a patient if you forget to reset your unit manually.

5. The top value on the response scale is 80. A top value of 80 was chosen because it gives you higher resolution in your results. If you reach a value of 80 in a test, the indicator light on the probe wand will start flashing. Once the light begins to flash, the testing sequence is complete.


6. False negative and false positive indications can be detected and avoided. A false negative indication occurs when a healthy tooth does not respond to maximum stimulus, i.e. a reading of 80. It can be caused by moisture on the surface of the tooth or by a metallic filling making contact with periodontal tissue. Other conditions that can cause a false negative reading are severe calcification of the pulpal tissue or other conditions that result in the degradation of the neural response.

A false positive indication occurs when the patient feels the stimulus even when the pulpal tissue is dead. This can be caused by a patient sensing the stimulation in an adjacent tooth which makes contact with the test tooth. The patient may also sense a periodontal neural stimulus rather than a stimulation of the pulpal nerve. The sensation of a periodontal response is quite different than a pulpal response, so we suggest that you compare the sensation of a tooth being pulp tested with the direct stimulation of periodontal tissue before you begin pulp testing.

7. The Vitality Scanner turns off automatically 10 to 15 seconds after contact is broken between the probe and the tooth. The automatic power-off feature saves battery life and permits you to install the unit at some distance from your chair since only the probe needs to be readily accessible.

8. The Vitality Scanner is powered by four 1.5V size AA alkaline batteries. The voltage is electronically stabilized so the function of the Vitality Scanner is independent of the battery voltage. The voltage stabilization eliminates the risk of faulty diagnosis due to the voltage fluctuations of depleted batteries. Two dots will begin to flash on the display when the batteries are getting low. You should change your batteries as soon as you see the two dots appear.

9. Accessories Model 2006 comes with a short and a long probe tip. Also included with each model is a grounding lead with lip clip. You may also purchase additional tips, grounding leads, and lip clips.

Using the Vitality Scanner

The following procedure is recommended for using the Vitality Scanner:

1. Before using the instrument for the first time, adjust the sweep rate to its lowest setting by turning the panel wheel to "1". Connect the grounding lead to the probe wand and insert the tip into the front of the probe wand.
2. Attach the lip clip to the patient's lip. Alternatively, the patient may choose to hold the lip clip in his/her hand.
3. Dry the tooth to be tested.
4. Dip the probe tip in a small amount of fluoride gel and place it on the side of the tooth. Do not make contact with the gum.
5. The light on the probe wand will turn on when you have achieved a good electrical contact. Hold the probe in contact with the tooth, and lift it from the tooth when the patient perceives the stimulus.
6. The response level is shown on the digital display.
7. After a delay of two seconds you can begin testing the next tooth. The display resets automatically when you make contact.
8. After you are finished testing your Vitality Scanner will turn off automatically after 10 to 15 seconds.

Normal response ranges are: Incisors 10–40
Bicuspids 20–50
Molars 30–70

Frequently asked questions about the Vitality Scanner

1. Why is the Vitality Scanner better than other units on the market?

The main reason why the Vitality Scanner is better is that it never jolts the patient like the other units do. Its patented contact sensing feature prevents the unit from turning on and increasing the output unless there is a good electrical contact with the tooth. With other units you never know if the current is actually flowing through the tooth. In some cases you may have mechanical contact with the tooth but not electrical contact. When this happens the output stimulus can get very high without the patient feeling anything. If contact is suddenly made at this high output level, the patient will receive a powerful and painful jolt.

Other advantages are the automatically increasing output which assures a gentle increase in the stimulus level and the automatic reset function,

which resets the output to its lowest level after lifting the tip away from the tooth. No longer do you have to remember to turn the output down before proceeding to the next tooth.

2. Can I use the unit while wearing gloves? Yes, a grounding lead with lip clip is delivered with each unit. The grounding lead provides a conductive path for the electrical current even if you are wearing gloves. One end clips onto the shaft of the probe wand, the other to an autoclavable stainless steel lip clip. This lip clip can be placed on the lip of the patient or can be held in the patient's hand.

3. What kind of conductive medium do you recommend between the tip and the tooth? We recommend topical fluoride gel. Toothpaste or an electrode paste may also be used, but the topical gel conducts electricity better than toothpaste and is generally more readily available to general practitioners than electrode paste. If you want to use toothpaste, gel type toothpastes work the best. However, please note that Sensodyne is a poor conductor and should be avoided.

4. Do I have to dry the tooth before using the Vitality Scanner? No, the Vitality Scanner is relatively insensitive to moisture. However, the tip should not contact the gingiva or a metallic filling which in turn could make contact with the gingiva or an adjacent tooth. All you need to do is place the electrode squarely on the side of the tooth.

5. How do I know when to change the batteries? The two decimal points for the display numbers with light up when the batteries run low.

6. What does the stimulus feel like from the patient's perspective? The sensation is a pulsating feeling which increases in strength. The patient will feel it long before it starts to hurt. Some people think it feels like a pressure, others like a warming or a tingling sensation. We recommend that you try it on yourself before using it on a patient so you have an idea of how the stimulus feels.

7. Can I test crowned teeth? Yes, we have a special mini-tip/crown test probe with a fine wire which can be used below the crown margin. The wire is insulated with only the tip surface exposed. It does not matter if the wire touches the gingiva. It works best, however, if the margin is open.

8. Can I use the device on a patient with a pacemaker? Although the current is very low we recommend that you do not use it on a patient with pacemaker.

9. Do I have to isolate the tooth being tested from adjacent teeth? Generally not. However, if there is a large metallic filling in contact with an adjacent tooth, isolation with a rubber dam may be necessary.

10. What rate of increase do you recommend? We recommend that you start with a low rate of increase until you get comfortable with the operation of the Vitality Scanner. A rate in the range of 5 to 7 is common for experienced practitioners.

11. Do I have to wait until the unit turns off before proceeding to the next tooth? No, you can proceed after waiting two seconds. The display will reset to zero as soon as you touch the next tooth with the tip.

Cleaning, Sterilization, and Maintenance

Sterilization

Probes and lip clips must be cleaned and sterilized before every use:

1. Clean with water and mild, non-abrasive, detergent such as dish washing liquid. Dry thoroughly.
2. Double wrap the probes and clips using a typical central service wrapping technique and place into autoclave unit avoiding contact with other instruments.
3. Operate sterilization cycle at 250°F for 30 minutes or 270°F for 25 minutes. Follow the instructions supplied with your autoclave unit for specific operating instructions.

Expected Life of Accessories

All accessories are reusable only for a limited number of procedures. This number is highly dependent on how and for what duration the instrument is used during each procedure and, therefore, will vary significantly from user to user. If the device function appears to be erratic, an accessory may have exceeded its useful life and should be replaced.

Cleaning

The Vitality Scanner 2006 enclosure, probe wand, cord, and grounding lead should be wiped down with a cloth soaked in any of the common non-abrasive, mild, cleaning liquids containing alcohol. Do not submerge the device or let liquid enter the enclosure.

Maintenance

To replace the batteries simply remove the two (2) screws and battery cover on the bottom of the unit. Replace with only high quality size AA alkaline batteries. During replacement, note proper battery alignment. Misalignment of battery cell polarity can lead to device malfunction. Replace the cover and screws. There are no other user serviceable items within the unit.

Technical Specifications

Dimension/Weight

Casing of unit: Weight: 2.50 lb. (1.13 kg)
Dimensions: Height: 2.0 inches/51 mm
Width: 4.3 inches/109 mm
Depth: 7.8 inches/198 mm

Classification: IEC 601-1 Internally Powered Equipment

Mode of Operation: IEC 601-1 Continuous Operation

Degree of safety of application in presence of a flammable anesthetic mixture:

Unit unsuitable for use in presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

Conditions for use: +10°C to +40°C
30% to 75% relative humidity

Storage and transport conditions: -10°C to +45°C
10% to 95% relative humidity

Origin: SybronEndo
1332 South Lone Hill Avenue, Glendora, CA 91740, USA

Model: Vitality Scanner 2006

Power Supply: 4 each, 1.5 volts AA alkaline batteries

Warranty Information

Returns and Warranty

The Vitality Scanner 2006 is designed for lasting durability and predictable results. The unit comes with a one-year warranty.

If the unit should fail to operate correctly, please contact your SybronEndo customer service representative at 1-800-346-3636. Outside the United States, please call customer service at 714-516-7979 or contact your local dealer or supplier to coordinate the returns/repairs with SybronEndo.

Warranty repairs can only be performed by SybronEndo or by authorized agencies with original factory parts. Any unauthorized repairs will void the warranty.

1. Use the original packaging when shipping or storing your Vitality Scanner.
2. Many problems can be solved over the phone. Do not hesitate to contact us if you experience difficulties when using the Vitality Scanner.
3. For returns, call us for a Return Good Authorization (RGA) number.
4. Damages incurred during shipping due to improper packaging are not covered by the warranty. If the original box and/or foam packing are not available, please call SybronEndo for packing instructions.
5. Mark all outer packaging with the RGA number, your name, address, and phone number.
6. Shipping costs are not covered by the warranty.

Disclaimer

The responsibility of SybronEndo, as a manufacturer of electrical/medical devices, extends to the technical safety features of the device only if maintenance and repairs are carried out by SybronEndo or by authorized agencies with original factory parts.

For safety reasons, this product should be used with accessories manufactured and sold by SybronEndo. Any use of non-authorized accessories or not following any of the instructions for use is done so at the operator's risk and voids the warranty.

SybronEndo does not assume any responsibility for incorrect diagnosis due to operator error or equipment malfunction.

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