VOCO

Alfacomp LC

GB

In accordance with DIN EN ISO 4049

CE 0482

Instructions for use

Product description:

Alfacomp LC is a light-curing, highly radiopaque hybrid composite, coloured grey to contrast more strongly with the natural dental hard tissue. Its physical stability and high abrasion resistance were developed to meet the special requirements of restorations bearing masticatory loads. **Alfacomp LC** contains 86 % w/w inorganic fillers (equivalent to 70.1 % by volume) in a methacrylate matrix (Bis-GMA, UDMA, TEGDMA) and is cured by blue halogen or LED light. It is available in syringes and Caps for direct application.

Shade:

Alfacomp LC is grey coloured.

Indications:

- Class I and II restorative treatment in the posterior region
- Core build-up

Type of application:

Preparation:

Clean teeth to be treated with a fluoride-free cleaning paste.

Moisture control:

Ensure an adequate moisture control. The use of a rubber dam is recommended.

Preparing the cavity:

The cavity should be prepared according to the rules for adhesive restorative treatment; preparation should be minimally invasive to protect healthy dental hard tissue. Then clean and dry the cavity. When applying bonding agents, observe any special characteristics (see relevant instructions for use).

Cavity lining:

Suitable protection for the pulp should be applied when working close to it and should be covered with a layer of stable cement if necessary.

Placing the matrix:

For cavities in the proximal region apply a matrix. It helps to use translucent matrices, wedged in the proximal region. Keeping separation to the minimum makes it easier to shape the proximal contacts and attach a matrix.

Bonding material:

Alfacomp LC is applied with a dentine enamel bond, using the adhesive technique. Any light-curing bonding materials can be used. Observe the relevant instructions for use during preparation and processing.

Processing:

Bring material to ambient temperature before application.

Application:

Apply and cure restorations exceeding 2 mm incrementally in layers.

Light-curing:

Conventional polymerisation devices are suitable for light-curing the material. Curing time is **10 seconds per 2 mm** increment when an LED/halogen lamp with a light output of at least 500 mW/cm² is used.

The light aperture of the light guide should be brought as close as possible to the restoration. Make sure that every area of the restoration is sufficiently exposed to the light. You are advised to expose the accessible proximal boxes to the light orally or vestibularly. If you are using metal matrices, remove them beforehand. Inadequate curing can cause discolouration and sensitivity.

Finishing:

You can finish and polish restorations using cooling, immediately after removing the shaping aids (e. g. fine or extra-fine diamond burs, abrasive discs). Finally, the edge of the restoration, or better still the whole tooth, should be fluoridated.

Warning/precautionary measures:

- **Alfacomp LC** contains Bis-GMA, UDMA, BHT, TEGDMA. **Alfacomp LC** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.
- Phenolic substances, such as products containing eugenol or thymol, result in curing failures in restorative composites. Consequently, you should avoid using zinc oxide eugenol cements or other materials containing eugenol in combination with restorative composites.
- Occlusally loaded restorations should be checked at least once a year so that any changes can be detected in good time.
- Each **Alfacomp LC** Caps unit should be used for one patient only and should be applied exerting uniform but not excessive pressure.

Storage:

Store at temperatures of 4 °C - 23 °C. If refrigerated, the material must be allowed to reach room temperature before use. To avoid exposure to light and potentially resulting polymerisation, syringes should be closed immediately after dispensing. Do not use after the expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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