Rebilda DC

Carefully read instructions prior to use

Rebilda DC is a dual-curing, radiopaque, flowable composite indicated for core build-ups.

Rebilda DC is available in the following shades:

- **blue** to visualize the transition between material and tooth substance, especially in molar and pre-molar areas;
- **dentin** shade for high esthetic demands, e.g. under full ceramic restorations or composite restorations with high translucency; and
- **white** as a possibility to visualize preparation margins and simultaneously meet esthetic demands.

Indications:

Core build-up of vital and non-vital teeth.

Contraindications:

Do not use in cases of known allergies to the following ingredients:
 - **Rebilda DC** contains Bis-GMA, UDMA, BHT and benzoyl peroxide

Directions:

A) Preparation

Take appropriate protective measures in areas close to the pulp.

B) Bonding material

Rebilda DC is used in the adhesive technique in combination with a light- or dual-curing adhesive which must be compatible with dual-curing composites. The product's instructions for use should be consulted to determine the correct preparation (etching technique) and procedure.

C1) Rebilda DC in the safety cartridge

Insert the **Rebilda DC** safety cartridge into the VOCO dispenser (type 3) or any other dispenser with a suitable piston. Due to technical reasons, slight differences in the filling level between the cartridge chambers can occur. In this case, remove cartridge cap and squeeze out the material until the material is delivered evenly from both chambers. Then, insert a mixing tip type 12 and lock with a 90° turn clockwise. If necessary, attach an additional intraoral tip type 3 on the top of the mixing tip. By pressing the lever, the material is extruded and automatically mixed in the correct ratio.

C2) Rebilda DC in the QuickMix syringe

Remove the cap from the syringe after turning it counter-clockwise. Attach a type 11 mixing tip. Match the nooks on the syringe and the mixing tip and fix the mixing tip with a 90° turn clockwise. If necessary, attach an application tip type 4 on top of the mixing tip. The material is mixed automatically in the correct ratio when dispensing and can be applied directly. All QM syringes come with a detachable piston. The piston can simply be pulled from the QM syringe and replaced without causing any material leakage. The plugs in the syringe chambers will stay in place, preventing any cross-contamination of the material in the syringe.

D) Apply the core build-up material.

Place a matrix around the prepared tooth and apply **Rebilda DC** directly from the safety cartridge/QuickMix syringe.

Light-curing

Transparent matrices are recommended. Place the tip of the polymerization lamp as close as possible to the surface of the core build-up and polymerize with halogen light for 40 sec. If the distance is more than 5 mm, the curing depth may be reduced. In case of a larger core build-up, apply in layers of 2 mm and cure each layer for 40 sec. After light-curing, the material can be finished immediately.

Chemical curing

The build-up can be finished 5 min after its application in the mouth without light-curing.

Rebilda DC has a working time of approx. 2 min. When using a prefabricated crown, apply **Rebilda DC** directly into the prefabricated crown, seat the crown on the prepared tooth, adjust its position and light-cure the margins for 40 sec.

After use, leave the mixing tip on the cartridge/QuickMix syringe or replace the cap.

Luting of fibre-reinforced resin posts (in combination with a dual-curing adhesive):

1. When luting root posts, apply the dual-curing adhesive according to its instructions for use.

2. **Using cartridge:**

When using intraoral tip type 5 as tip on top of the mixing tip, apply **Rebilda DC** directly into the root canal. For this purpose the intraoral tip can be pre-bent by hand, if necessary. If not using intraoral tip type 5, proceed as follows: After preparing the post according to the manufacturer's instructions, wet post with **Rebilda DC**. Do not apply **Rebilda DC** directly into the root canal or use a lentulo! Insert the post making sure to leave a slight excess.

Using QuickMix syringe:

It is recommended to apply **Rebilda DC** directly into the root canal using thin intraoral tips type 1. Do not use a lentulo! After preparing the root pin according to the manufacturer's instructions, it can be additionally moistened with **Rebilda DC**. Insert the post making sure to leave a slight excess.

3. To fix the post and cure the dual-curing adhesive applied to the occlusal surface, light-cure for at least 40 sec.
4. For core build-up, follow the instructions given in D) Application of the core build-up material.

Caution:

- For complete polymerization, avoid any contact of the adhesive and **Rebilda DC** with phenolic substances (e.g. eugenol).
 - Do not use bonding agents which are not compatible with chemically curing or dual-curing composites.
 - Insulate the core build-up before making temporary crowns/bridges (good moistening is often sufficient) to avoid bonding to these materials.
 - Avoid contact with the oral mucosa.
- The physical values of **Rebilda DC** correspond to EN 24049/ISO 4049.

Storage:

Close **Rebilda DC** cartridges/QuickMix syringes after use. Store in a cool area. Do not expose **Rebilda DC** to direct sunlight or the chairlight. Store at temperatures between 39 °F - 73 °F = 4 °C - 23 °C. Do not use materials after expiration date. Do not sanitize the mixing tips in ultrasonic sterilizers.

Product Ordering Information:

Intro Kit dentin: 50gm cartridge, 8ml Futurabond DC , Mixing tips type 12, extruder gun type 3, accessories	1395
Pack 50gm cartridge Rebilda DC dentin , mixing tips type 12	1396
Pack 50gm cartridge Rebilda DC blue , mixing tips type 12	1397
Pack 50gm cartridge Rebilda DC white , mixing tips type 12	1398
Intro Kit dentin: QuickMix syringe, 8ml Futurabond DC , mixing tips type 11, accessories	1402
Pack 10gm QuickMix syringe Rebilda DC dentin , mixing tips type 11	1403
Pack 10gm QuickMix syringe Rebilda DC blue , mixing tips type 11	1404
Pack 10gm QuickMix syringe Rebilda DC white , mixing tips type 11	1405

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. VOCO recognizes its responsibility to replace products if proven to be defective. VOCO does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith. Descriptions and data constitute no warranty of attributes and are not binding.

CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

For questions or comments, please call 1-888-658-2584.

**Keep this material out of reach of children.
For dental use only.**

An explanation of the symbols used in labeling can be found at www.voco.dental/symbols

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VOCO GmbH Phone +49 (4721) 719-0